

# Creator™ DNA Cloning Kits User Manual



Cat. Nos. 631614, 631615  
PT3460-1 (PR7Y2438)  
Published 28 November 2007

## Table of Contents

---

<b>I. Introduction</b>	<b>4</b>
<b>II. List of Components</b>	<b>12</b>
<b>III. Additional Materials Required</b>	<b>13</b>
<b>IV. Creator™ DNA Cloning System</b>	<b>14</b>
A. Before you start	14
B. Creator DNA Cloning Procedure	14
C. Colony PCR	16
<b>V. Typical Results</b>	<b>17</b>
<b>VI. Troubleshooting Guide</b>	<b>18</b>
<b>VII. References</b>	<b>21</b>
<b>Appendix A: Creator™ Donor &amp; Control Vector Maps</b>	<b>22</b>
<b>Appendix B: Creator™ Donor Vector MCSs</b>	<b>25</b>
<b>Appendix C: Competent Cells</b>	<b>27</b>

---

**Table of Contents** *continued*

---

**List of Figures**

Figure 1. <i>LoxP</i> sequence.	4
Figure 2. Flow chart of the Creator DNA Cloning System.	5
Figure 3. Flow chart of tagging using pDNR-Dual.	8
Figure 4. Typical Plating Result	16
Figure 5. The Creator system easily generates 10 different constructs in one day.	17
Figure 6. Amplification across a recombination juncture.	18
Figure 7. Typical test results for successful recombination.	20
Figure 8. Map of pDNR-1r Donor Vector.	22
Figure 9. Map of pDNR-1r-Luc Control Vector.	22
Figure 10. Map of pDNR-Dual Donor Vector.	23
Figure 11. Map of pDNR-Dual-Luc Control Vector.	23
Figure 12. Map of pDNR-CMV Donor Vector.	24
Figure 13. Map of pDNR-CMV-LacZ Control Vector.	24
Figure 14. MCS of pDNR-1r Donor Vector.	25
Figure 15. MCS of pDNR-Dual Donor Vector.	25
Figure 16. MCS of pDNR-CMV Donor Vector.	26

**List of Tables**

Table I. Creator Acceptor Vectors	9
Table II. Creator Donor Vectors	10
Table III. Recommended Competent Cells	27


## I. Introduction

The **Creator™ DNA Cloning Kits** provide a revolutionary system for transferring a target gene directly into multiple expression vectors, without the need for time-consuming subcloning steps. With the Creator DNA Cloning Kits, you can study novel protein interactions, tetracycline-regulated expression, signal transduction, retroviral expression, fluorescent protein tagging, and many other applications, simultaneously.

### Cre-*lox* site-specific recombination

The Creator System uses Cre-*loxP* site-specific recombination to catalyze the transfer of a target gene from a Donor Vector plasmid containing your gene of interest to an Acceptor Vector, a plasmid containing regulatory elements of the desired host expression system (Figure 2). Cre, a 38-kDa recombinase protein from bacteriophage P1, mediates recombination between or within DNA sequences at specific locations called *loxP* sites (Sauer, 1994; Abremski *et al.*, 1984). These sites consist of two 13-bp inverted repeats separated by an 8-bp spacer region that provides directionality to the recombination reaction (Figure 1). The 8-bp spacer region in the *loxP* site has a defined orientation which forces your gene to be transferred in a fixed orientation and reading frame.

Donor Vectors contain two *loxP* sites, which flank the 5' end of the MCS and the 5' end of the open reading frame for the chloramphenicol resistance gene (Cm<sup>r</sup>; Figure 2). Donor Vectors also contain the ampicillin gene for propagation and selection in *E. coli*, and the sucrose gene from *B. subtilis* (*SacB*) for selection of correct recombinants. Acceptor Vectors contain a single *loxP* site, followed by a bacterial promoter, which drives expression of the chloramphenicol marker after Cre-*lox*-mediated recombination. The gene of interest, once transferred, will become linked to the specific expression elements for which the Acceptor Vector was designed. Furthermore, if the coding sequence for the gene of interest is in frame with the upstream *loxP* site in the Donor Vector, it will automatically be in frame with all peptides in the Acceptor Vector. Therefore, you only need to determine the correct reading frame once, and your target gene will always be transferred in the



**Figure 1.** *LoxP* sequence, showing reading frame.

# I. Introduction *continued*




Figure 2. Flow chart of the Creator™ DNA Cloning System showing a representative Donor Vector, pDNR-1r.

## I. Introduction *continued*

---

correct reading frame and orientation in the Acceptor Vector after recombination and selection.

### **Creator™ DNA cloning & expression system**

Figure 2 provides an overview of the Creator DNA Cloning & Expression System. To transfer your gene from the Donor Vector into any Acceptor Vector, simply incubate the Donor Vector containing your gene of interest and an Acceptor Vector with Cre Recombinase (see Tables I & II for a list of Creator Acceptor and Donor Vectors). Cre binds to the *loxP* sites on both the Donor Vector and Acceptor Vector, cleaves the DNA, and covalently attaches itself to the DNA. Then Cre catalyzes strand exchange and ligation of the DNA so that the gene is transferred from the Donor Vector into the Acceptor Vector. As a result, a recombinant construct is created that expresses your gene of interest in the desired host system. Chloramphenicol and sucrose selection lets you harvest desired recombinant colonies that contain a directionally correct gene insert. Clones containing the remaining Donor Vector, without your gene insert, will express *SacB*, and therefore, cannot be grown on media containing sucrose.

### **Transfer target genes easily**

Traditional cloning practices require several days of tedious restriction enzyme digestion, fragment purification, and re-ligation procedures. The Creator System is so simple that in just one day you can create multiple constructs that are ready for immediate use in functional studies.

In separate tubes, combine appropriate Acceptor Vectors with a Donor Vector, containing your gene of interest, and Cre Recombinase. Incubate tube(s) at room temperature for 15 minutes. Next, the enzyme is heat-killed for 5 minutes at 70°C. Transform competent cells with an aliquot of the Creator reaction mixture. After 30–60 minutes in SOC or LB medium, plate cells on agar plates containing chloramphenicol and sucrose for selection of desired recombinants. After colonies are selected, you can immediately prepare DNA for further studies. If desired, recombinant plasmids can be further propagated in either chloramphenicol or the antibiotic that is appropriate for the resistance marker carried by the Acceptor Vector.

### **Creator gives you access to multiple expression systems**

With Creator you can place your gene within multiple systems for functional analysis. Clontech offers Acceptor Vectors for many of our most powerful systems (Table I). The Matchmaker™ Acceptor Vectors enable you to discover novel protein-protein interactions using the Matchmaker Two-Hybrid System 3. The Living Colors™ Acceptor Vectors allow you to generate C-terminal fusions of your target gene to a fluorescent tag for protein localization studies. For dose-dependent inducible expression studies of your protein, transfer your gene into an Acceptor Vector that is compatible with both

---

## I. Introduction *continued*

---

Tet and the retroviral RevTet™ Gene Expression Systems. The IRES bicistronic expression Acceptor Vectors have an IRES (internal ribosomal entry site) sequence and a constitutive CMV promoter to produce a bicistronic message for high expression of your protein in mammalian cells. The pLP-PROtet-6xHN Acceptor Vector (Cat. No. 631201) yields high inducible protein levels in bacterial cells for easy purification with TALON® resin. For more information on these and other Creator vectors, please refer to Tables I and II or visit the Clontech page at [www.clontech.com](http://www.clontech.com).

Complete Systems compatible with Creator are available for several of our gene expression systems. These complete kits offer a streamlined approach to expressing your target gene. See Section VIII, Related Products for a current list.

### **Creator Cloning Kits**

The Creator™ pDNR Cloning Kit (Cat. No. 631615) is a reformulated version of our original Creator Cloning Kit. The Donor Vector in this kit, pDNR-1r is designed to increase recovery of recombinant vectors containing insert. The pDNR Cloning Kit can be used to transfer a target gene into any of our wide range of expression Acceptor Vector in the presence of Cre recombinase.


### **Donor Reporter Vectors generate reporter constructs in 15 minutes**

Our Donor Reporter Vectors, pDNR-SEAP, and pDNR-LacZ allow you to rapidly generate reporter constructs for any gene expression system. Using the standard 15 minute recombination reaction, Cre recombinase mediates the transfer of the reporter from the Donor Vector to any expression Acceptor Vector. The resulting reporter Acceptor Vectors are suitable for many uses, including transfection efficiency controls or quantifiers of induction of gene expression. These reporters are also ideal for normalizing your data.

### **Ready-made Creator™ Libraries Ñor generate your own**

The Creator™ SMART™ cDNA Libraries were developed to serve as the foundation for the Mammalian Gene Collection (MGC) project, a joint effort of the National Institutes of Health (NIH) and the National Cancer Institute (NCI). This project aims to facilitate the search for and isolation of your target gene by providing researchers with a full set of inexpensive, full-length clones and sequences from human and other mammalian sources. Using a full-length sequencing pipeline, the project produces 5,000 highly accurate full-length sequences per year. These sequences are analyzed by a bioinformatics group, and then made available to the research community through GenBank and the MGC website (<http://mgc.nci.nih.gov>). By providing the biomedical research community with inexpensive, full-length clones in a rapid-cloning format, the MGC project hopes to accelerate the functional analysis of genes identified by mammalian genome projects. Currently, nearly 1,000 full-length, sequenced-verified cDNAs that are cloned into

## I. Introduction *continued*



**Figure 3. Flowchart of tagging using pDNR-Dual.** The pDNR-Dual Donor Vector contains a splice donor (SD) site directly downstream of the Multiple Cloning Site (MCS). The SD site, which is transferred from pDNR-Dual along with the gene of interest, mediates the fusion of the gene to the tag in the Acceptor Vector through intron splicing. As a result, a transcript is created that expresses the tag as a 3' fusion to your gene of interest in a eukaryotic host system.



## I. Introduction *continued*

**TABLE I. CREATOR™ ACCEPTOR VECTORS**

<b>Acceptor Vector</b>	<b>Promoter/Features</b>	<b>Functional Application</b>
<b>pLP-GADT7</b>	ADH1/GAL4 activation domain	Express fusions to GAL4 AD to study protein interactions through two-hybrid screening.
<b>pLP-GBKT7</b>	ADH1/GAL4 DNA-binding domain	Express fusions to GAL4 DNA-binding domain to study protein interactions through two-hybrid screening.
<b>pLP-AcGFP1-C</b>	CMV/ fusions of gene to terminal of AcGFP	Express fusions to AcGFP ( <i>Aequora C-coerulescens</i> green fluorescent protein) to study localization of the protein of interest in live cells; no dyes or cofactors required.
<b>pLP-IRESneo</b>	CMV/IRES, neo selection marker	Constitutive mammalian expression with single transcript for both gene of interest and neo selection marker.
<b>pLP-TRE2</b>	Inducible tet-responsive expression vector	High-level, regulated promoter for expression in mammalian cells.
<b>pLP-RevTRE</b>	Inducible tet-responsive promoter in retroviral expression vector with hyg selection	High-level, regulated retroviral expression in mammalian cells.
<b>pLP-LNCX</b>	CMV promoter in retroviral expression vector with neo selection	Constitutive retroviral expression.
<b>pLP-CMVneo</b>	CMV promoter with neo selection	Constitutive mammalian expression
<b>pLP-CMV-Myc</b>	CMV/C-terminal fusions to Myc	Express fusions to Myc tag for detection
<b>pLP-CMV-HA</b>	CMV/C-terminal fusions to HA	Express fusions to HA tag for detection
<b>pLP-PROTet-6xHN</b>	Inducible tet-responsive expression vector	High-level, regulated promoter for expression in bacteria through two-hybrid screening
<b>pLP-BacPAK9</b>	Prokaryotic promoter baculoviral construct	Constitutive baculoviral expression

\* Only compatible with genes cloned in pDNR-Dual Donor Vector.

## I. Introduction *continued*

**TABLE II. CREATOR DONOR VECTORS**

Donor Vector	Promoter/Features	Functional Application
pDNR-1r	T7 promoter, M13 F/R primer sites, SacB selection	T7 RNA polymerase primer/promoter site upstream of MCS for <i>in vitro</i> transcription/translation of gene of interest
pDNR-Dual	T7 promoter, M13 forward site, SD site, 6xHN c-ter tag, SacB selection	T7 RNA polymerase primer/promoter site upstream of MCS for <i>in vitro</i> transcription/translation of gene of interest C-ter tagging by intron splicing in Eukaryotes and has built in 6xHN tag at C-ter for bacteria.
pDNR-CMV	CMV promoter, M13 F/R primer sites, T7 promoter, SacB selection	T7 RNA polymerase primer/promoter site upstream of MCS for <i>in vitro</i> transcription/translation of gene of interest, CMV promoter for expression testing in mammalian cells prior to transfer
pDNR-LIB	T7 promoter, M13 F/R primer sites, SacB selection	smaller pDNR 1 vector, T7 RNA polymerase primer/promoter site upstream of MCS for <i>in vitro</i> transcription/translation of gene of interest, designed for library construction using SMART™ Library Construction Kit
pDNR- LacZ pDNR- SEAP	SacB selection	Donor Vector with Reporter genes for enzymatic quantitation or cell staining assay

the pDNR-LIB Vector from the Creator SMART cDNA Libraries are available through the I.M.A.G.E. Consortium. The libraries can be quickly screened, and your full-length cDNA clones isolated by standard procedures. The target gene is then efficiently transferred into a variety of expression vectors for functional analysis. The Creator Libraries are available from numerous cancer and normal tissues from various sources.

The Creator™ SMART™ cDNA Library Construction Kit (Cat. No. 634903) provides a dependable method for producing high-quality, cDNA libraries compatible with Creator Systems specific for your research needs. SMART technology makes it possible to generate full-length, directionally cloned cDNA Libraries from nanograms of total or poly A<sup>+</sup> RNA. Using the Creator Cloning System, isolated clones from finished libraries can be transferred directly to expression Acceptor Vectors for functional analysis without the need for subcloning.

### Advantages over other Cre-*lox* systems

In the Creator System, two *loxP* sites in the Donor Vector flank the MCS and chloramphenicol open reading frame, so only this small region of the Donor Vector is transferred to the Acceptor Vector. In other systems, there is

---

## I. Introduction *continued*

---

only one *loxP* site in the Donor Vector causing both the donor and Acceptor Vectors to fuse into one large plasmid. Thus, single *loxP* site-based Donor Vectors are not compatible with retroviral and IRES Acceptor Vectors due to size constraints. The small expression vector produced by the Creator System is easier to use and is ideal for a variety of downstream applications, including retroviral expression.

Creator does not require PCR cloning methods so there is no need to sequence the entire gene insert. Your sequence remains intact, an especially important feature for genes undergoing functional studies. Other systems require you to sequence the full insert because their procedures use PCR-based cloning, which may introduce errors due to the low fidelity of some thermostable polymerases. Another advantage of the Creator System is that you can transform any standard competent cell lines, such as DH5 $\alpha$  cells, after performing Cre-mediated recombination. In addition, Cre Recombinase is the only enzyme necessary for the recombination reaction, and it does not require additional cofactors, nor does it have a preference for linear over supercoiled DNA.

The Creator DNA Cloning Kits include Cre Recombinase, a Donor Vector, a Control Donor Vector, 10X BSA, and 10X Cre Reaction Buffer. Each kit includes sufficient reagents to perform 10 reactions. In addition, Cre Recombinase can also be purchased separately (Cat. No. 631614; sufficient for 20 reactions). Clontech offers a variety of competent cells and a wide range of Acceptor Vectors.

## II. List of Components

---

Store vectors at  $-20^{\circ}\text{C}$ . Store all other components at  $-70^{\circ}\text{C}$ .

After the first use, store Cre Recombinase at  $-20^{\circ}\text{C}$  in a non-frost-free freezer.

Creator DNA Cloning Kits each include sufficient reagents to perform 10 reactions as described in this User Manual.

### Creator pDNR Cloning Kit (Cat. No. 631615)

- 20  $\mu\text{g}$  pDNR-1r Vector (500 ng/ $\mu\text{l}$ )
- 20  $\mu\text{g}$  pDNR-1r-Luc Control Vector (500 ng/ $\mu\text{l}$ )
- 10  $\mu\text{l}$  Cre Recombinase
- 100  $\mu\text{l}$  10X Cre Reaction Buffer
- 100  $\mu\text{l}$  10X BSA (1 mg/ml)
- pDNR-1r Vector Information Packet (PT3616-5)

### Cre Recombinase (Cat. No. 631614)

- 20  $\mu\text{l}$  Cre Recombinase
- 100  $\mu\text{l}$  10X Cre Reaction Buffer
- 100  $\mu\text{l}$  10X BSA (1 mg/ml)

### III. Additional Materials Required

---

The following materials are required but not supplied:

- **Acceptor Vectors**

(For ordering information, see [www.clontech.com](http://www.clontech.com).)

- **Electrocompetent or chemically-competent cells**

See Table III. (For ordering information, see [www.clontech.com](http://www.clontech.com).)

**Note:** "One shot" TOP10 cells are not compatible with the Creator Cloning System.

- **Chloramphenicol stock solution (Cm)**

(30 mg/ml in 100% ethanol; 1000X.) Store at  $-20^{\circ}\text{C}$ .

- **Ampicillin stock solution (Amp)**

(100 mg/ml in  $\text{H}_2\text{O}$ ; 1000X.) Store at  $-20^{\circ}\text{C}$ .

- **SOC Medium**

20 g/L Bacto-tryptone

5 g/L Bacto-yeast extract

0.5 g/L NaCl

Add 10 ml of a 250 mM solution of KCl, and adjust the pH to 7.0. Autoclave. Prior to use, add 5 ml of 2 M  $\text{MgCl}_2$  (autoclaved).

- **LB broth**

10 g/L Bacto-tryptone

5 g/L Bacto-yeast extract

5 g/L NaCl

Adjust pH to 7.0 with 5 N NaOH. Autoclave.

- **LB/amp agar plates**

Prepare LB as described above. Prior to autoclaving, add agar (15 g/L). After autoclaving, cool to  $55^{\circ}\text{C}$  and add ampicillin (100  $\mu\text{g}/\text{ml}$ ; final concentration). Pour into 10-cm plates and let harden. Then invert plates and store at  $4^{\circ}\text{C}$ .

- **LB/Cm/sucrose agar plates**

Prepare LB as described above. Prior to adjusting pH, add sucrose (7% w/v). Add agar (15 g/L) and autoclave for 20 minutes maximum (longer autoclaving times may cause sucrose to burn). After autoclaving, cool to  $55^{\circ}\text{C}$  and add chloramphenicol (30  $\mu\text{g}/\text{ml}$ ). Pour into 10-cm plates, and let harden. Then invert plates and store at  $4^{\circ}\text{C}$ .

## IV. Creator™ DNA Cloning Protocol

### A. Before you start

- Although the reading frame of the insertion is only significant when transferring into Acceptor Vectors that express fusion proteins, we strongly encourage you to clone your gene of interest in frame with the Donor Vector's upstream *loxP* site to ensure compatibility with all of our Acceptor Vectors. For detailed Donor Vector maps and multiple cloning site diagrams, see Appendices A & B.
- When using the pDNR-Dual Cloning Kit for the addition of 3' tags to a gene of interest, the sequence cloned into pDNR-Dual must be in frame with the SD site, and lack stop codons and a 3' UTR for correct expression of the protein-tag fusion.
- Cre recombinase requires heat to efficiently inactivate the enzyme following the reaction. The enzyme retains activity even when reactions are stored at 4°C.
- Sucrose and chloramphenicol selection is only required for initial selection of correct recombinants. After selection, recombinant clones can be grown in the same medium used for propagation of the Acceptor Vector.
- For best recombination efficiency, we recommend using Nucleo-Spin® Plasmid Miniprep or Maxiprep Kits to obtain highly pure plasmid DNA (see [www.clontech.com](http://www.clontech.com) for ordering information).
- Verify plasmid concentration by gel electrophoresis before performing the recombination reaction. Please note that spectrophotometer readings can be inaccurate in particular when RNA is column purified.
- Please see Appendix C for a list of recommended competent cells.
- We strongly recommend that you perform a control reaction using the Control Donor Vector provided in each kit. Performing a control reaction will verify that your system is working properly.

### B. Creator™ DNA Cloning Procedure

1. Subclone your gene of interest into the Donor Vector using standard methods. (For more information, consult Sambrook *et al.*, 2001.)
2. Prepare the Creator reaction mixture as follows:

200 ng	<b>Donor Vector</b>
200 ng	<b>Acceptor Vector</b>
2 µl	<b>10X Cre Reaction Buffer</b>
2 µl	<b>10X BSA (1 mg/ml)</b>
1 µl	<b>Cre Recombinase</b>

Add deionized H<sub>2</sub>O to bring volume up to 20 µl.

---


## IV. Creator™ DNA Cloning Protocol *continued*

---

**Notes:**

- We recommend performing a parallel control reaction with the appropriate Control Vector as described above (Section A).
  - The first time you use Cre Recombinase, thaw rapidly from  $-70^{\circ}\text{C}$  by holding the tube between your fingers until almost all of the ice is melted. Then place the tube on ice.
  - Once thawed, the Cre Recombinase should be stored at  $-20^{\circ}\text{C}$  in a non-frost-free freezer. The Cre solution will remain a liquid and will not need to be thawed prior to use.
3. Mix well by gently tapping the tube. Spin briefly.
  4. Incubate at room temperature ( $22\text{--}25^{\circ}\text{C}$ ) for 15 min.  
**Note:** Do not extend incubation past 15 min; competing recombination reactions, which do not generate desired recombinants, can reduce the yield of your desired recombinants.
  5. Stop the reaction by heating tube at  $70^{\circ}\text{C}$  for 5 min.
  6. Transform competent cells with 1  $\mu\text{l}$  of reaction mixture. See Appendix C for a list of recommended competent cells and volumes to use.  
**Note:** Competent cells should give  $>1 \times 10^8$  cfu/ $\mu\text{g}$ . If not, replace with a fresh sample of cells.
  7. Grow cells at  $37^{\circ}\text{C}$  for 60 min in SOC medium (or LB).
  8. Plate 100  $\mu\text{l}$  of transformation on a 10-cm LB-agar plate containing 30  $\mu\text{g}/\text{ml}$  chloramphenicol, and 7% sucrose (w/v).
  9. Centrifuge at 6,000 rpm for 1 min.
  10. Aspirate off 800  $\mu\text{l}$ .
  11. Resuspend in 100  $\mu\text{l}$  of SOC medium (or LB) and plate on a 10-cm LB agar plate containing 30  $\mu\text{g}/\text{ml}$  chloramphenicol, and 7% sucrose (w/v). Incubate overnight.  
**Note:** It is important to allow the inoculum to soak into the plate thoroughly before incubating overnight.
  12. The next day, the plate should contain a mixture of larger colonies and smaller colonies. Pick larger colonies for screening by colony PCR (Figure 4). Smaller colonies typically contain a mixture of Donor and Acceptor Vector DNA. Analyze clones by restriction digest to confirm your final construct. Plasmid(s) can be further propagated in either chloramphenicol or the antibiotic that is appropriate for the resistance marker of the Acceptor Vector.

## IV. Creator™ DNA Cloning Protocol *continued*



**Figure 4. Typical plating result.** Arrows indicate examples of larger colonies that are good choices for selection and screening by colony PCR.


### C. Colony PCR

1. Create a 0.4  $\mu$ M primer master mix. 50  $\mu$ l of this mix will be required per colony screened. For recombinations involving the acceptor vectors pLP-CMV-HA or pLP-CMV-Myc, use a combination of primers PCP-1 and PCP-L. For all other acceptor vectors use PCP-1 and PCP-2. Primer sequences are:  
PCP-1: GCTCACCGTCTTTCATTGCC  
PCP-2: TCCGCTCATGAGACAATAACC  
PCP-L: TGTATCTTATCATGTCTGGATC
2. Aliquot 50  $\mu$ l of primer master mix to wells of a microtiter plate or strip. Pick large colonies from the plate in step 12 and pipet up and down several times to resuspend the bacteria. Use 25  $\mu$ l of this master mix of bacteria and primers to resuspend the contents of one well of a Sprint Advantage 96 plate (Cat. No. 639550). Alternatively, other methods for "colony PCR" may be employed. It is important however, to add broth (eg. LB or SOC) containing 30  $\mu$ g/ml Chloramphenicol or another antibiotic that is appropriate for your plasmid to the resuspended colony within ten minutes. This will ensure the viability of antibiotic-resistant bacteria.
3. Perform the PCR and analyze as described in the troubleshooting guide (No or few colonies > Recombination reaction failed: Steps 3–5).



## V. Typical Results

Figure 5 shows typical results using the Creator DNA Cloning System to generate multiple constructs, each containing a target gene. In separate tubes, pDNR-Luc Control Vector was combined with Cre Recombinase and ten different Acceptor Vectors. After a brief recombination, DH5 $\alpha$  competent cells were transformed with an aliquot of reaction mixture and recombinants were selected. With only three colonies selected from each transformation, we obtained the desired recombinants containing luciferase from all 10 reactions. Out of 30 colonies picked in total, all but two contained the proper insert. With the Creator DNA Cloning Kits, 95–100% of analyzed clones will contain your desired recombinant construct.



**Figure 5. The Creator™ system easily generates many constructs in one day.** Luciferase was transferred into each of ten Acceptor Vectors. Three clones for each construct were then screened by restriction digest to confirm transfer of the luciferase gene. All but two of the 30 colonies contained the correct recombinant construct. P = parental vector. M = marker. \*Discontinued products as of 3/3/2005. Data shown was obtained prior to 3/3/2005, pLP-AcGFP1-C is now available as a replacement for pLP-EGFP-C1, pLP-ECFP-C1, and pLP-EYFP-C1.

## VI. Troubleshooting Guide

The simplicity of the Creator DNA Cloning System makes its use fairly straightforward. If you do not achieve typical results, this guide may help you to determine the source of the problem.

### No or few colonies

Colonies of varying size

- Sometimes distinct large and small colonies may be seen. In such cases, we generally find it best to pick the larger colonies. Small colonies tend to be mixtures of recombinant plasmid and Acceptor Vector.

Low transformation efficiency

- Check transformation efficiency. You should obtain  $>1 \times 10^8$  cfu/ $\mu$ g; otherwise use new, fresh competent cells.

Too much Cre used in the reaction


- To ensure that an accurate volume of Cre recombinase is added to the reaction mix, use a 2- or 10- $\mu$ l pipettor.

Cell type used

- While in general there is no specific requirement to use a given *E. coli* strain, there is some variability in overall efficiencies in different cell types. Thus, if you have low colony numbers or get colonies that are mixtures (see above), then trying a different cell type can help. Test both electro- and chemically competent versions of different cell types (e.g., DH10B, XL1-Blue, Fusion-Blue, or DH5 $\alpha$ ) for best results.

Incubation of reaction exceeded 15 minutes

- **Do not incubate Creator Reaction Mixture longer than 15 minutes.** Competing recombination events can reduce the yield of your desired recombinants.
- Ensure that Cre recombinase is heat-killed immediately following incubation. We have found that the reaction is relatively temperature insensitive, and can even occur on ice. For this reason, it is critical that you closely monitor the reaction time and proceed to the next



**Figure 6. Amplification across a recombination junction.** The two primers, PCP-1 and PCP-2, prime from the chloramphenicol marker and prokaryotic promoter, respectively. Since these two elements are only linked through recombination, successful PCR amplification indicates a successful recombination reaction.

## VI. Troubleshooting Guide *continued*

step in the protocol immediately, rather than let the reaction sit on ice.

### Poor plasmid DNA quality

- For best results, we recommend using NucleoSpin Plasmid Products to obtain highly pure, miniprep plasmid DNA. For ordering information, see [www.clontech.com](http://www.clontech.com).
- To accurately determine plasmid DNA concentration, analyze an aliquot of your miniprep, alongside a Molecular weight standard, on a 0.8–1.0% agarose/EtBr gel. Photograph the gel and compare the band intensity of the plasmid to the band intensities of the mass standards to quantify. Please note that spectrophotometer readings can be inaccurate in particular when DNA is column-purified.

### Recombination reaction failed

1. To test the success of the recombination reaction, set up a PCR reaction mixture using TITANIUM™ Taq (Cat. No. 639208) and primers (PCP-1 and PCP-2) that amplify across a recombination juncture (see Figures 2 and 6). Please note: Use Primer PCP-L instead of PCP-2 if using pLP-CMV-Myc or pLP-CMV-HA as Acceptor Vectors.
2. Set up the PCR reaction in a PCR tube as follows:

Test Sample	Negative Control	
18.5 µl	19.5 µl	PCR-Grade Water
2.5 µl	2.5 µl	10XTITANIUM Taq PCR Buffer
1 µl		Creator Reaction Mixture
1 µl	1 µl	PCP-1 Primer (10 µM)
1 µl	1 µl	PCP-2 Primer (10 µM)
0.5 µl	0.5 µl	50X dNTP Mix (10 mM ea. of dATP, dCTP, dGTP, dTTP)
0.5 µl	0.5 µl	50XTITANIUM Taq DNA Polymerase
<b>25 µl</b>	<b>25 µl</b>	<b>Total volume</b>

Primer PCP-1: 5'-GCTCACCGTCTTTCATTGCC-3'

Primer PCP-2: 5'-TCCGCTCATGAGACAATAACC-3'

Primer PCP-L: 5'-TGTATCTTATCATGTCTGGATC-3'


3. Commence thermal cycling using the following parameters. These are general guidelines for use with hot-lid thermal cyclers; the optimal parameters may vary with different thermal cyclers.
  - 94°C for 2 min
  - 25 cycles:
 

94°C	15 sec
58°C	30 sec
72°C	30 sec
  - 72°C for 5 min

## VI. Troubleshooting Guide *continued*

4. Transfer a 5- $\mu$ l sample of your PCR reaction to a fresh tube and add 1  $\mu$ l of 5X stop/loading buffer. Analyze your sample(s) by electrophoresis on a 2.0% agarose/ethidium bromide gel, along with suitable DNA size markers.
5. **Expected results:** The reaction should produce a major fragment of 356 bp. No bands should be generated in the negative (i.e., no DNA template) control.

**Note:** This test is also a very simple way to screen colonies for recombinants.



**Figure 7.** Typical test results for successful recombination.

---

## VII. References

---

Abremski, K. & Hoess, R. (1984) Bacteriophage P1 site-specific recombination. Purification and properties of the Cre recombinase. *J. Biol. Chem.* **259**:1509–1514.

Creator DNA Cloning & Expression System (April 2000) *Clontechniques XI*(2):7–11.

Sambrook, J., Fritsch, E. F. & Maniatis, T. (2001). *Molecular Cloning: A Laboratory Manual*, Cold Spring Harbor Laboratory Press (Cold Spring Harbor, NY).

Sauer, B. (1994) Site-specific recombination: developments and applications. *Curr. Opin. Biotechnol.* **5**:521–527.


### For additional reading

Barabino S. M. L. & Keller, W. (1999) Last but not least: regulated Poly(A) Tail Formation. *Cell* **99**:9–11.


Guo, F., Gopaul, D. N. & Van Duyne, G. D. (1997) Structure of Cre recombinase complexed with DNA in a site-specific recombination synapse. *Nature* **389**:40–46.

Lilley, D. M. J. (1997) Site-specific recombination caught in the act. *Chem. Biol.* **4**:717–720.

## Appendix A: Creator™ Donor and Control Vector Maps




**Figure 8. Map of pDNR-1r Donor Vector.** The Donor Vector contains two *loxP* sites, which flank the 5' end of the MCS and the 5' end of the chloramphenicol open reading frame (*Cm<sup>r</sup>*). The Donor Vector also contains the sucrose gene from *B. subtilis* (*SacB*) for selection of correct recombinants, and an ampicillin resistance gene for propagation and selection in *E. coli*. Each Acceptor Vector contains a *loxP* site, followed by a bacterial promoter, which drives expression of the chloramphenicol marker after recombination. The gene of interest, once transferred, will become linked to the specific expression elements for which the Acceptor Vector was designed. In addition, if the gene of interest is cloned in frame with the upstream *loxP* site in the Donor Vector, it will automatically be in frame with all peptides in the Acceptor Vector following recombination. Sequence and digest information is available, and can be downloaded from our web site at [www.clontech.com/clontech/techinfo/vectors/](http://www.clontech.com/clontech/techinfo/vectors/).




**Figure 9. Map of pDNR-1r-Luc Control Vector.** This vector is similar to the Donor Vector, above, except the luciferase gene has been cloned into the MCS.

## Appendix A: Creator™ Vector Maps *continued*




**Figure 10. Map of pDNR-Dual Donor Vector.** pDNR-Dual contains a splice donor (SD) site directly downstream of the Multiple Cloning Site (MCS). When combined with a specialized Acceptor Vector containing a splice acceptor (SA) site, a recombinant expression construct is generated containing an artificial intron (consisting of the chloramphenicol marker and one *loxP* site), which is spliced out by the eukaryotic host's transcriptional machinery. As a result, a transcript is created that expresses the tag as a 3' fusion to your gene of interest. Sequence and digest information is available, and can be downloaded from our web site at [www.clontech.com/clontech/techinfo/vectors/](http://www.clontech.com/clontech/techinfo/vectors/).




**Figure 11. Map of pDNR-Dual-Luc Control Vector.** This vector is similar to the Donor Vector, above, with the exception that the luciferase open reading frame has been cloned into the MCS.

## Appendix A: Creator™ Vector Maps *continued*




**Figure 12. Map of pDNR-CMV Donor Vector.** pDNR-CMV contains the immediate early CMV promoter for expression testing in mammalian cells prior to gene transfer, two *loxP* sites (flanking the 5' end of the MCS) and the 5' end of the chloramphenicol open reading frame (*Cm<sup>r</sup>*). The Donor Vector also contains the sucrase gene from *B. subtilis* (*SacB*) for selection of correct recombinants, and an ampicillin resistance gene for propagation and selection in *E. coli*. Each Acceptor Vector contains a *loxP* site, followed by a bacterial promoter, which drives expression of the chloramphenicol marker after recombination. Sequence and digest information is available, and can be downloaded from our web site at [www.clontech.com/clontech/techinfo/vectors/](http://www.clontech.com/clontech/techinfo/vectors/).




**Figure 13. Map of pDNR-CMV-LacZ Control Vector.** This vector is similar to the pDNR-CMV Donor Vector, shown above, with the exception that the *lacZ* gene has been cloned into the MCS.



## Appendix B: Creator™ Donor Vector MCSs




**Figure 14. MCS of pDNR-1r Donor Vector.** Unique restriction sites are shown in bold. The MCS is shown in frame with the *loxP* site. The last four nucleotide bases of the *loxP* site can be seen at the left hand side of the map in bold. If the coding sequence for the gene of interest is in frame with the upstream *loxP* site in the Donor Vector, it will automatically be in frame with all peptides in the Acceptor Vector.



**Figure 15. MCS of pDNR-Dual Donor Vector.** Unique restriction sites are shown in bold. The MCS is shown in frame with the *loxP* site (frame 1). The last four nucleotide bases of the *loxP* site can be seen at the left hand side of the map in bold. If the coding sequence for the gene of interest is in frame with the upstream *loxP* site in the Donor Vector, it will automatically be in frame with any 5' peptides in the Acceptor Vector. If the coding sequence for the gene of interest is in frame with the SD site in the Donor Vector, it will automatically be in frame with any 3' tags in the Acceptor Vector. The *Apa* I site in the MCS is methylated and will only cut in a dam<sup>-</sup> host strain.

## Appendix B: Creator™ Donor Vector MCSs *continued*



**Figure 16. MCS of pDNR-CMV Donor Vector.** Unique restriction sites are shown in bold. The MCS is shown in frame with the *loxP* site (frame 1). The last four nucleotide bases of the *loxP* site can be seen at the left hand side of the map in bold. If the coding sequence for the gene of interest is in frame with the upstream *loxP* site in the Donor Vector, it will automatically be in frame with any 5' peptides in the Acceptor Vector.

## Appendix C: Competent Cells

Creator cloning results may vary depending on the competent cells used. At Clontech, we have used the Creator System successfully with many commercially available competent cells. In our experience, the cells listed in Table I have yielded the best results.

**TABLE III. RECOMMENDED COMPETENT CELLS**

<b>Competent Cells</b>	<b>Volume of Cells (μl)<sup>a</sup></b>
<b>Chemically Competent</b>	
Max Efficiency® DH5α (Invitrogen Corporation Cat. No. 18297-010)	50
Fusion-Blue™ Competent Cells (Cat. No. 636700)	50
NovaBlue Singles™ (Novagen No. 70181-3)	50
XL1-Blue (Stratagene No. 200236)	100
<b>Electrocompetent</b>	
Supercharge EZ10 Electrocompetent Cells (Cat. No. 636756)	40
ElectroMAX™ DH10B™ (Invitrogen Corporation Cat No. 18290-015)	20

<sup>a</sup> Per 1 μl of recombination reaction mix.

## Notes

---

### Notice to Purchaser

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without written approval of Clontech Laboratories, Inc.

The **Creator™** Technology is based on the process of *in vitro* Cre-LoxP recombination. Clontech is the assignee of U.S. Patent Nos. 6,410,317; 6,977,165; 6,838,285 and other patents pending covering Creator vectors and the selection process as it relates to the production of recombinant clones. Clontech has chosen not to exercise its right to impose license fees on the in-house use of the Creator Technology. However a royalty-bearing license is required on contract services and sale or distribution of clones made in the Creator System format. For information on using Creator Technology for commercial purposes, please contact a licensing representative by phone at 650.919.7320 or by e-mail at [licensing@clontech.com](mailto:licensing@clontech.com).

### For **Living Colors®** Products

AcGFP1, DsRed, HcRed, AsRed, AmCyan, ZsGreen, ZsYellow and their variants:

These products are the subjects of pending U.S. and foreign patents.

**Not-For-Profit Entities:** Orders may be placed in the normal manner by contacting your local representative or Clontech Customer Service at 650.919.7300. At its discretion, Clontech grants Not-For-Profit Entities a non-exclusive, royalty-free, personal, limited license to use this product for non-commercial life science research use only. Such license specifically excludes the right to sell or otherwise transfer this product, its components or derivatives thereof to third parties. No modifications to the protein coding sequence may be made without express written permission from Clontech. Any other use of this product requires a license from Clontech. For license information, please contact a licensing representative by phone at 650.919.7320 or by e-mail at [licensing@clontech.com](mailto:licensing@clontech.com).

**For-Profit Entities** wishing to use this product are required to obtain a license from Clontech. For license information, please contact a licensing representative by phone at 650.919.7320 or by e-mail at [licensing@clontech.com](mailto:licensing@clontech.com).

The CMV promoter is covered under U.S. Patent Nos. 5,168,062, and 5,385,839 assigned to the University of Iowa Research Foundation.

Practice of the **two-hybrid system** is covered by U.S. Patent Nos. 5,283,173, 5,468,614, and 5,667,973 assigned to the Research Foundation of the State University of New York. Purchase of any Clontech two-hybrid reagent does not imply or convey a license to practice the two-hybrid system covered by these patents. Commercial entities purchasing these reagents must obtain a license from the Research Foundation of the State University of New York before using them. Clontech is required by its licensing agreement to submit a report of all purchasers of two-hybrid reagents to SUNY Stony Brook. Please contact the Office of Technology Licensing & Industry Relations at SUNY Stony Brook for license information (Tel: 631.632.9009; Fax: 631.632.1505).

Use of the **IRES sequence** is covered by U.S. Patent No. 4,937,190 and is limited to use solely for research purposes. Any other use of the IRES sequence requires a license from Wisconsin Alumni Research Foundation.

Use of the **Tetracycline controllable expression systems** (the "Tet Technology") is covered by a series of patents including U.S. patent Nos. 5,464,758 and 5,814,618, which are proprietary to TET Systems Holding GmbH & Co. KG. Academic research institutions are granted an automatic license with the purchase of this product to use the Tet Technology only for internal, academic research purposes, which license specifically excludes the right to sell, or otherwise transfer, the Tet Technology or its component parts to third parties. Notwithstanding the above, academic and not-for profit research institutions who's research using the Tet Technology is sponsored by for profit organizations, which shall receive ownership to all data and results stemming from the sponsored research, shall need a commercial license agreement from IP Merchandisers in order to use the Tet Technology. In accepting this license, all users acknowledge that the Tet Technology is experimental in nature. TET Systems Holding GmbH & Co. KG makes no warranties, express or implied or of any kind, and hereby disclaims any warranties, representations, or guarantees of any kind as to the Tet Technology, patents, or products. All others are invited to request a license from TET Systems Holding GmbH & Co. KG prior to purchasing these reagents or using them for any purpose. Clontech is required by its licensing agreement to submit a report of all purchasers of the Tet-controllable expression system to IP Merchandisers, Inc. For license information, please contact:

Hans Peter Kneubuehl  
TET Systems Holding GmbH & Co. KG  
Im Neuenheimer Feld 582  
69120 Heidelberg  
Germany

Tel +49 6221 588 04 00

Fax +49 6221 588 04 04

eMail: [kneubuehl@tet-systems.de](mailto:kneubuehl@tet-systems.de)

or use our electronic licensing request form via [http://www.tetsystems.com/main\\_inquiry.htm](http://www.tetsystems.com/main_inquiry.htm)

This product is also sold under patent sublicense FOR RESEARCH PURPOSES ONLY. Licenses for commercial manufacture or use under this patent may be obtained directly from Harvard University.

DH10B™, DH5 $\alpha$ ™, ElectroMAX™ are trademarks and Max Efficiency® is a registered trademark of Invitrogen Corporation.

NucleoSpin® is a registered trademark of MACHEREY-NAGEL GmbH and Co., KG.

Singles™ is a trademark and Novagen® is a registered trademark of Novagen, Inc.

Clontech, the Clontech Logo and all other trademarks are property of Clontech Laboratories, Inc. unless noted otherwise. Clontech is a Takara Bio Company. ©2007 Clontech Laboratories, Inc.